





## **FT100B**

## **INFRARED DIGITAL THERMOMETER\***

- LCD display wtih 3½ figures
- Temperature unit of measurement: Celsius degrees (°C) and Fahrenheit (°F)
- Measurement range: 32.0 42.9 °C (89,6 109.2 °F)
- Saves last 10 performed measurements
- Powered by not rechargeable battery, PC 2032 type (1pc)
- The batteries last about 2500 uses
- Size: 105x35x35mmNet weight: 36gr
- · Acoustic alarm: beep sound for ON-OFF
- Lifetime: 3 years
- Complete with protective plastic case

DIR/QM - Rev. 02\_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.