



## BOR81145N

### ATLS 2.0 RESCUE RUCKSACK\*

- ATLS rescue rucksack
- Made of tear-resistant and waterproof nylon
- Features reflective stripe on padded shoulder straps for improved visibility
- The transparent front window is suitable for displaying of a logo
- The large front pocket allows you to safely store even large medical devices
- The backpack is equipped with soft rear pads to reduce back fatigue
- The shoulder straps are equipped with a sternum closure for easier transport
- The backpack is equipped with a rigid shockproof polypropylene tray with central oxygen cylinder holder support
- The 4 internal coloured pouches have a transparent front and are removable
- Dimensions: about 32x25x53cm
- Volume: about 35 litres
- Weight: about 3.3kg
- Maximum gross weight: 30kg
- The ATLS backpack is sold empty but can be set up on specific customer request

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.