



HS-201Q1

SPHYGMOMANOMETER PROFESS. I

- Professional 1 tube sphygmomanometer
- Big manometer with scale up to 300 mm/Hg (70 mm diameter)
- Cuff with Velcro closing
- Deflation valve with screw regulator or rapid deflation button
- Turning manometer that allows the use to both left- and right-handed users
- Practical protection and transport case
- Is a medical device CE
- The device can be used by trained personnel only or doctors

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.