



LD-588

AUTOMATIC/DIGITAL SPHYGMOMANOMETER*

- 52x59 mm LCD display
- Up to 60 measurements storage capacity
- Oscillometric measurement method
- Pressure measurement scale: 40-260 mm/Hg
- Pulse measurement scale: 40-160 beats/min
- Accuracy: ± 3 mm/Hg (pressure) $\pm 5\%$ (pulse)
- Automatic pressurization
- Quick depressurization with valve
- Power supply: 4 x 1.5V AA batteries (not included)
- Dimensions: 127x99x50 mm
- Weight: approx. 255 g without batteries
- It is a CE medical device
- Measurements are to be carried out by qualified healthcare professionals

DIR/QM – Rev. 02_04/09/2024

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.