



x50 pz

207010

GISAFE SURGICAL MASK BLACK 50 PCS

- Filtering masks 98% with 3 layers - TYPE IIR with elastic bands.
- Lightweight and comfortable surgical masks suitable for hospitals, dental and medical clinics.
- Adaptable nose clip
- Bacterial filtration efficiency (BFE) = 98%.
- Made of non-woven, latex-free and glass fibre-free material
- Colour: black
- Dimensions: approx. 17,5 x 9,5 cm
- Pack of 50 pcs
- Class I medical device in accordance with MDR 2017/745

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.