



**SC125003**

## **SUCTION UNIT CATHETER CH12 WHITE**

- Straight-tipped bronchial suction probe, single-use and sterile.
- Blunt tip, atraumatic, with two elliptical staggered side holes.
- Made of non-toxic medical-grade PVC, soft, with satanized surface for easy insertion and removal. Phthalate-free.
- Storage: in a cool, dry place, away from heat sources, in the original sales packaging. Temperature limits: 10-45° C.
- Dimensions for Ch. 12. Length: 50 cm. Outer diameter: 4 mm. Color: white.
- Medical device in accordance with the Medical Device Directive 93/42/EEC, class of risk Is (I sterile)

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.