



FS-128

RESCUE ISOTHERMIC BLANKET*

- Protection and rescue isothermic blanket
- Prevents patients' hypothermia whether traumatized, burned or post-operated
- Double function: warming and cooling
- Silver side, if positioned toward the body, isolates and reflects, preserving about 80% of body warmth. Thus protects patient from hypothermia, coldness, freezing and shock
- Golden side: placed toward patient body and silver one, reflects sun rays maintaining body freshness
- Non sterile
- Dimensions: about 160x210 cm

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.