



1333122

SILICONE RESUSCITATOR PEDIATRIC SIZE

- Transparent silicone manual resuscitator, autoclavable up to 121°C
- Volume 500 ml
- Complete with 360° rotating connection patient's valve
- Safety valve for overpressure regulation
- Additional inlet for dynamometer
- Peep valve inlet
- Supplied with mask no. 3, reservoir and oxygen tube
- CE certified in class IIa

User manual available on the [eIFU](http://eifu.com) webpage

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.