



**BOR87220R**

## **RUCKSACK MASTER TROLLEY EMPTY**

- Trolley Rucksack Master empty, suitable for carrying materials in between hospital departments
- Made of red tearproof fabric
- 2 refractive strips, pull-out trolley handle, pull-out padded shoulder-straps
- 2 frontal pockets which might be equipped with additional material, 4 removable coloured pouches with transparent front
- Sold empty (additional material on request)
- Dimensions: about 31x37x58 (h) cm
- Weight: about 5,5 kg

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.