









## **BOR87040**

## **OXIBAG BLUE\***

- Padded oxygen cylinder case
- Waterproof and tear-proof high-strength PVC fabric
- The shape and padding allow the oxygen cylinder to be safely contained, protected and transported
- Can hold 1 or 2 l oxygen cylinders
- Dimensions: approx. 54x23x24 cm
- Weight: approx. 1.1 kg

DIR/QM - Rev. 01\_26/07/2023

For further information and/or technical data related to the product, please refer to the Oscar Boscarol Company.