



**SC185003**

## **SONDA ASPIRAZIONE 18 CH ROSSO**

- Bronchial suction catheter with straight tip, single-use and sterile.
- Blunt, atraumatic tip with two elliptical staggered side holes.
- Made of non-toxic medical grade PVC, soft and with a satinised surface for easy insertion and removal.
- Phthalate-free.
- Dimensions:
  - Length: 50 cm.
  - Size: CH 18.
  - Outer diameter: 6.0 mm.
- Colour: red.
- EC Medical Device

*DIR/QM – Rev. 02\_27/07/2023*

Per ulteriori informazioni e/o dati tecnici relativi al prodotto, fare riferimento alla Oscar Boscarol srl